



A FORUM FOR SUPPORTERS OF THE DISTRIBUTED NETWORK PROTOCOL

## Distributed Network Protocol

This certificate confirms that the device described below has successfully complied with the requirements as outlined by the DNP 3.0 Technical Committee's DNP 3.0 Conformance Testing procedure. This certificate entitles the bearer to promote the below stated device as having fully passed the DNP 3.0 Subset Definitions v2.6 rev1.

**Manufacturer:** Siemens


**Device Model:** SIPROTEC 4, 7SC80

**Firmware Number:** V01.00.05

**Level Successfully Passed:** Level 2

**Testing House:** The testing house is KEMA, with the DNP logo also present.

**Date:** November 23, 2012



Tester:

A handwritten signature in blue ink, appearing to read 'P.H.S. Ermens'.

P.H.S. Ermens